INSTRUCTIONS FOR USE

Ambu® Neuline® Cup Neurology surface electrodes

Single-use

Abridged version - en

1. Important information – Read before use
These safety instructions carefully follow using the Ambu® Neuline™ Cup Neurology surface electrodes. The instruction for use may be updated without further notice. Copies of the current instructions are available upon request. Please be aware that this document does not explain or discuss clinical procedures. They describe only the basic operation and precautions related to the operation of the electrodes.

There is no warranty on the Ambu Neuline Cup Neurology surface electrodes.

1.1. Intended use
Ambu Neuline Cup electrodes are intended for use with recording and monitoring equipment for the purpose of recording biopotential signals for example: Electroencephalography (EEG) and nerve potential signals.

1.2. Indications for use
Ambu Neuline Cup electrodes are indicated when there is a need for neurophysiological examinations of the brain to investigate abnormalities in a patient.

1.3. Intended patient population
Ambu Neuline Cup electrodes can be used for adults (age 22 and up*) and for pediatric (from both through the age of 19).

*Age definitions are based on FDA Guidance: Leveraging Existing Clinical Data for Extrapolation to Pediatric Uses of Medical Devices (Sep 19, 2016).

1.4. Intended user/ use environment
Ambu Neuline cup electrodes are intended for used by in consultation with nurse or physician.

The electrode is for use in a hospital/clinic setting.

1.5. Contra indications


Note: **CAUTIONS**

1.6. Clinical benefits
Facilitates neurophysiological examination through detecting and transmitting the signal to the recording and monitoring equipment.

1.7. Warnings and cautions

- WARNINGS**
- Surface electrodes should be applied only to intact, clean skin (e.g., not over open wounds, lesions, infected, or inflamed areas).
 - For single use only. Reuse of the product may cause contamination, leading to infections.
 - Do not use the electrode for electrical stimulation as this may lead to tissue discomfort/skin burn.
 - Do not use the electrode on patient with body weight less than 3.5 kg as this may lead to skin irritation.

- CAUTIONS**
- To minimize the risk of skin irritation, when replacing electrodes find a new site of skin on which to place them.
 - Electrodes should only be used by or in consultation with a health care provider familiar with their proper placement and use.
 - Do not soak, rinse, or sterilize the electrodes as these procedures may leave harmful residues or cause malfunction of the electrodes.
 - During surgical procedures electrodes should be placed as far as possible from the electro-surgical unit to avoid electrical coupling.
 - US Federal law restricts this device to sale by or on the order of a licensed health care practitioner.

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1.8. Undesirable side effects
Skin irritation, skin redness, skin injuries such as skin burn and skin abrasion.

1.9. General notes
If you are using this device as a result of its use, a serious incident may occur. Please report it to the manufacturer and to your national authority.

2. Device description
Ambu Neuline Cup electrodes are single-use surface electrodes consisting of cup, crimp lock followed by a shrink tube and a pre-attached lead wire with a connector which enables connection to a neurological acquisition system.

Ambu Neuline Cup electrodes are available with pre-attached lead wire of 100 cm (40”), 150 cm (60”) and 200 cm (80”) lengths, and with a touch proof 1.5-mm safety connector. Pre-attached lead wire is available in 12 colours.

The packaging configurations of the Ambu Neuline Cup electrode includes 10, 15, 23, 25 and 27 electrodes per primary packaging.

The packaging configurations of the Ambu Neuline Cup electrode includes 10, 15, 23, 25 and 27 electrodes per primary packaging.

Symbol indication	Description	Symbol indication	Description
	Medical Device		Date of Manufacture Country of Manufacture
	Electrical Safety Type BF Applied Part		Prescription use only
	UK Conformity Assessed		UK Responsible Person
	Importer (For products imported into Great Britain only)		

A full list of symbol explanations can be found on ambu.com/symbol-explanation.

4. Product use
• The correct placement of electrodes is typically determined by the hospital/physician and/or local protocol/procedures.
• Prepare the scalp using an abrasive. Apply gel, paste or adhesive on the cup electrode. Attach the cup to the scalp.
• Connect the lead wire to the monitoring device.
• After monitoring, remove the electrode from the skin.
• The electrodes are disposable and must be disposed of according to local hospital procedures.

5. Technical specifications
5.1. Specifications
Duration of use
Duration of use is up to 24 hours.

Storage temperature
The product shall be stored at 10 to 25 °C (50 to 77° F).

- Correct placement of electrodes is typically determined by the hospital/physician and/or local protocol/procedures.
- Prepare the scalp using an abrasive. Apply gel, paste or adhesive on the cup electrode. Attach the cup to the scalp.
- Connect the lead wire to the monitoring device.
- After monitoring, remove the electrode from the skin.
- The electrodes are disposable and must be disposed of according to local hospital procedures.

- Correct placement of electrodes is typically determined by the hospital/physician and/or local protocol/procedures.
- Prepare the scalp using an abrasive. Apply gel, paste or adhesive on the cup electrode. Attach the cup to the scalp.
- Connect the lead wire to the monitoring device.
- After monitoring, remove the electrode from the skin.
- The electrodes are disposable and must be disposed of according to local hospital procedures.

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