


CARDIOLOGY ACCESSORIES

Features and Benefits

- Fixes electrode leadwires on patients during Holter and Stress test
- Reduces movement artefacts on the ECG
- Minimises the risk of electrodes falling off
- MR safe
- Adhesive size: Ø 55 mm
- PVC free
- Latex free
- Biocompatible



Skin Fix

Ref: 2703/50

Units/pouch	Units/inner box	Units/outer box
50	500	3000

Features and Benefits

- Allows skin preparation before electrode placement
- Gently abrasive to reduce skin impedance
- Adhesive on one side, to use with one finger
- Small size fits finger tip
- Size: 19 x 28 mm
- PVC free



Skin Prep

Ref: 2121M

Units/pouch	Units/inner box	Units/outer box
500	1500	15000

CARDIOLOGY ACCESSORIES

Features and Benefits

- Connects 3 mm pin cables to tab electrodes
- Ensures secure connection
- Easy to use
- Reusable
- PVC free



Snap Clip

Ref: 0334M

Units/pouch	Units/inner box	Units/outer box
10	N/A	1000

Features and Benefits

- Connects 4 mm banana pin to snap or tab electrodes
- Ensures secure connection
- Allows quick and easy one-hand application
- Reusable
- PVC free



Multi Clip

Ref: 2952/10

Units/pouch	Units/inner box	Units/outer box
10	N/A	N/A

Ambu

Ambu Australia Pty Ltd
Unit 5, 2 Daydream Street
Warriewood
NSW 2102
Tel: 1300 233 118
aussales@ambu.com
www.ambuaustralia.com.au

AA22 0095 V01
Always read the Instructions for Use. These products are only available for sale to healthcare professionals.

