

Ambu[®] Man Basic



Key features

- Cost efficient Basic Life Support (BLS) training manikin
- Robust and reliable
- Realistic training of chest compression
- Visualisation of the compression depth
- Ambu hygienic system for mouth-to-mouth and mouth-to-nose ventilation
- Can be used for bag-mask ventilation
- Placement of Guedel tubes (when using the Guedel face piece)
- Realistic ventilation resistance
- Visualisation of the correct ventilation volume
- Visible chest rise during ventilation
- Anatomically correct closed Torso
- Easy AED training

AmbuMan Basic

The AmbuMan Basic is a cost efficient training manikin that is designed to fulfill all requirements for the BLS algorithm training.

The closed shaped torso has anatomically correct landmarks. Trainees can also learn and train the use of Automatic External Defibrillators (AED) and the placement of pads.

The AmbuMan Basic enables an optimal technique of chest compression and ventilation with realistic resistance and expiration through mouth and nose.

The unique hygienic system prevents cross contamination.

The correct compression depth and required ventilation volume can be monitored by the feedback indicator on the manikin.

Realistic ventilation

The trainees will feel realistic airway resistance and the chest will start rising when the head is correctly tilted and air blown into the manikin.

Look, listen and feel according to BLS guidelines.

Cost efficient and reliable

The Ambu Man Basic is made to be cost efficient and reliable for BLS training.

Class training

For easy transportation and set up the AmbuMan Basic comes in a soft carry bag and packed as either a single manikin or in a set of four manikins.

Specifications

Weight Torso	6 kg
Length	70 cm
Width	37 cm
Weight 4 x Torso and Bag	30 kg
Max. Ventilation volume	1200 ml
Max. compression depth	7.5 cm

Materials

Basic unit	ABS (Acrylonitrile butadiene styrene)
Head	Polyvinyl Chloride, hard
Chest skin	Polyvinyl Chloride, soft
Face piece	Polyvinyl Chloride, soft
Headbag, disposable	Polyethylene
Training suit	Polyester
Carry bag	PVC-coated Nylon

Order information

Item no	Description
A268 005 000	AmbuMan Basic inc. training jacket, carry case, 1 face piece and 25 head bags
A268 006 000	AmbuMan Basic 4 pack inkl. 4 AmbuMan Basic Torso with training jacket, carry case, 4 face pieces and 100 head bags
Accessories	
234 000 702	Ambu head bags, disposable (set of 100)
234 000 703	Ambu face piece (Set of 5)
234 000 725	Ambu Guedel face piece



AmbuMan Basic



Hygienic System



Compression indicator

Ambu Australia Pty Ltd
Unit 5, 2 Daydream St
Warriewood, NSW 2102
www.ambuaustralia.com.au

Tel: 1300 233 118
aussales@ambu.com
AA21 0010 V01