

# Ambu<sup>®</sup> Man School



## Key features

- Ambu Hygienic system to prevent cross contamination
- Compact and portable training manikin
- easy training of chest compressions
- Closed torso
- Easy placement of AED pads
- 10 manikins for easy class training

## AmbuMan School

Ambu Man School is a training manikin which offers human appearance dedicated to the training of all procedures for basic life support resuscitation

The manikin can be used for realistic training of external chest compressions and the placement of AED pads in order to get used in the handling of AED's.

The unique hygienic system prevents cross contamination when demonstrating and training of the principals of ventilation through mouth and nose.

The AmbuMan School intended to be used for education of pupils in schools under the supervision of an instructor/ teacher.

### Guideline conform

The AmbuMan School allows to perform CPR training following the latest international CPR.

### Compact and portable

For an easy transportation & set up possibility the Ambu Man School has a low weight and compact packaging.

### Easy AED training

Due to the closed shaped torso, defibrillation pads can be placed and simulated AED training can be performed realistically.

## Specifications

|                        |        |
|------------------------|--------|
| Weight                 | 0,8 kg |
| Length                 | 45 cm  |
| Max. compression depth | 6 cm   |

## Material

|            |                         |
|------------|-------------------------|
| Straps     | Nylon                   |
| Base plate | Polycarbonate           |
| Bellow     | Thermoplastic Elastomer |
| Chest      | Polypropylene           |
| Head       | Silicone                |
| Face piece | Polyvinylchloride       |
| Lung bag   | Polyvinylchloride       |
| head bag   | Polyethylene            |

## Order information

| Item no.     | Description                            |
|--------------|--|
| A268 007 000 | AmbuMan School (Set of 10)             |
| 268 000 571  | Face piece AmbuMan School (Set of 10)  |
| 268 000 572  | Head bag AmbuMan School (Set of 100)   |
| 268 000 573  | Accessories AmbuMan School (Set of 10) |



Closed torso



Placement of AED pads



Ambu hygienic system



Accessories

**Ambu Australia Pty Ltd**  
 Unit 2, 1 Prosperity Parade  
 Warriewood, NSW 2102  
 Tel. 1300 233 118 Fax 1300 156 526  
 aussales@ambu.com  
 www.ambuaustralia.com.au