

DATASHEET

The choice for nasal endoscopy, flexible laryngoscopy and swallowing assessments.

The single-use aScope™ 4 RhinoLaryngo Slim video endoscope is designed with a small 3.0 mm outer diameter, to help minimise patient discomfort during procedures. It delivers high quality imaging with no reprocessing or repairs required and is easily available.

HD imaging helps you perform procedures confidently

The combination of the high image quality from the aScope™ 4 RhinoLaryngo Slim with the full-HD Ambu® aView™ 2 Advance monitor delivers excellent imaging. You can count on clear sharp imaging, which makes it easier to identify anatomical structures - and it allows sound recording when connected to the aView™ 2 Advance Generation 2.

Single-use efficiency streamlines your workflow

Whether you're on a consult or in the clinic, time is always of the essence. aScope™ 4 RhinoLaryngo is always available and easy to take with you streamlining your entire workflow.

Ambu® aScope™ 4 RhinoLaryngo Slim



KEY BENEFITS

High quality imaging on the full-HD monitor
Ambu® aView™ 2 Advance

It's an always available solution that supports increasing
patient volumes

Improves workflow without the hassle of cleaning and
repair

Sound and slow motion recording is available when
connected to the aView™ 2 Advance Generation 2.

A brand-new endoscope with the same high quality
imaging at every procedure

SPECIFICATIONS

Optical system

Field of view	85°
Direction of view	0° (forward viewing)
Depth of field	6 - 50 mm
Illumination method	LED

Insertion portion

Bending section	130° up, 130° down
Insertion cord diameter	3.0 mm (0.12")
Distal tip diameter	3.5 mm (0.14")
Maximum diameter of insertion portion	3.5 mm (0.14")
Working length	300 mm (11.8")

Operating environment, storage & transportation

Temperature - transportation	10 - 40° C (50 to 104° F)
Temperature - operating environment	10 - 40° C (50 to 104° F)
Temperature - storage	10 - 25° C (50 to 77° F)
Relative humidity	30 - 85%
Atmospheric pressure	80 - 109 kPa

Sterilisation

Method of sterilisation	ETO
-------------------------	-----

ORDERING SPECIFICATIONS

Product Code	Quantity
510001000 - aScope 4 RhinoLaryngo Slim	5



Ambu® aScope™ 4 RhinoLaryngo Slim and Ambu® aView™ 2 Advance



aScope™ 4 RhinoLaryngo Slim is comfortable to hold, thanks to its light-weight and ergonomic design

Ambu FOREVER FORWARD

Ambu Australia Pty Ltd | Ambu New Zealand Pty Ltd

Unit 5, 2 Daydream Street, Warriewood, NSW 2102

Tel (AUS): 1300 233 118 | Tel (NZ): 0800 804 079

Email: aussales@ambu.com

www.ambuaustralia.com.au

AA23 0133 V02. Always read the Instructions for Use.

These products are only available for sale to healthcare professionals.

CE 2797